

Solid Phase Extraction APPLICATIONS MANUAL

TABLE OF CONTENTS

Solvents, Solvents Mixtures and Use of Non-chlorinated Solvents	3
Reagents	4
How to Prepare Solutions and Buffers	5-6
Extraction Hints.....	7

Drug Testing Applications

Amphetamines In Urine	8
Anabolic Steroids in Urine	9-10
Barbiturates In Urine	11
Barbiturates In Urine (Reduced Solvent Volume).....	12
Basic Drugs for HPLC Analysis.....	13
Benzodiazepines In Urine	14-15
Benzodiazepines In Serum or Plasma	16
Benzodiazepines (CLEAN-UP®)	17
Beta Blockers In Urine	18
Beta Agonists In Urine	19
Carboxy THC In Urine	20
Carboxy THC In Urine (Reduced Solvent Volume)	21
Carisoprodol and Meprobamate For GC Or GC/GC/MS	22
Clonazepam & 7-Aminoclonazepam In Urine	23
Cocaine and Benzoyllecgonine In Urine	24
Cocaine and Benzoyllecgonine In Serum, Plasma Or Whole Blood For GC Or GC/MS	25
Cocaine and Benzoyllecgonine In Serum, Plasma Or Whole Blood For HPLC	26
Cocaine and Its Metabolites From Meconium For GC Or GC/MS	27
DHEA, Testosterone and Epitestosterone In Urine	28
Fentanyl and Analogues In Urine	29
Flunitrazepam and Metabolites In Urine	30-31
Gabapentin In Serum, Plasma Or Whole Blood	32
Gamma-Hydroxybutyrate (GHB) In Urine Without Conversion To Gamma-Butyrolactone (GBL)	33
Gamma-Hydroxybutyrate (GHB) In Blood Procedure	34
Gamma-Hydroxybutyrate (GHB) In Blood, Urine, Vitreous Or Tissue Without Conversion To Gamma-Butyrolactone (GHL)	35-36

Drug Testing Applications (Continue)

Glycopyrrolate (ROBINUL) From Equine Urine	37
Ketamine in Urine	38
Lysergic Acid Diethylamide (LSD) In Serum, Plasma Or Whole Blood	39
Lysergic Acid Diethylamide (LSD) In Urine	40
Manual Method For Immunoassay	41
Methaqualone In Urine	42
Methadone In Urine	43
Methylmalonic Acid in Serum Or Plasma	44
Nicotine And Cotinine In Urine Or Serum	45
Opiates In Urine - Oxime TMS Procedure For GC Or GC/MS	46-47
Opiates In Human Urine - Propyl Derivatives For GC Or GC/MS	48-49
Opiates (Free) In Serum, Plasma Or Whole Blood	50
Phencyclidine In Urine	51
Propoxyphene In Urine For GC Or GC/MS	52
Psilocin In Urine For GC Or GC/MS	53-54
Sertraline And Desmethylsertraline In Serum, Plasma Or Whole Blood	55
THC And Carboxy THC In Whole Blood For GC/MS	56-57
Therapeutic And Abused Drugs (Acid/Neutral/Basic) In Urine, Serum, Plasma Or Whole Blood	58-60
CLEAN SCREEN® DAU Forensic Applications	60
Tricyclic Antidepressants In Serum and Plasma For HPLC	61-62

Environmental Applications

4'4-Methylenedianiline in Serum	63
Chlorophenoxy Acid Herbicides in Water	64
Pesticides in Water	65
Polychlorinated Biphenyls in Pond Water	66
Polynuclear Aromatic Hydrocarbons in Pond Water	67

Miscellaneous

Abused Drugs in Canine or Equine Urine	68
Extraction of Tear Gas	69

SOLVENTS

Acetone: HPLC Grade
 Acetonitrile (CH₃CN): HPLC Grade
 Chloroform (CHCl₃): HPLC Grade
 Distilled or Deionized Water (DI H₂O): 5 < pH < 7
 Ethyl Acetate (EA): HPLC Grade
 Hexane: HPLC Grade
 Isopropyl Alcohol (IPA): HPLC Grade
 Methanol (CH₃OH): HPLC Grade
 Methylene Chloride (CH₂Cl₂): HPLC Grade

SOLVENT MIXTURES

Acetone/Hexane (1/99)
 Acetonitrile/DI H₂O (20/80)
 Ethyl Acetate/IPA (75/25)
 Hexane/Ethyl Acetate (50/50), (75/25)
 Methanol/DI H₂O (80/20)
 Methanol/DI H₂O (70/30)
 Methanol/DI H₂O (10/90)

USE OF NON-CHLORINATED ELUTION SOLVENTS

In response to environmental concerns over the use of chlorinated compounds in the laboratory, UCT offers these suggested non-chlorinated elution solvents. The recommended parameters have been used successfully on Worldwide Monitoring[®] columns by our customers across the country and may be routinely used as an alternative to chlorinated elution solvents. You may however see subtle differences on certain compounds due to solubility effects.

Assay	Chlorinated	Non-chlorinated
Opiates Propoxyphene	CH ₂ Cl ₂ /IPA/NH ₄ OH (78/20/2)	EA/IPA/NH ₄ OH (90/6/4)
Cocaine/BE Amphetamines	CH ₂ Cl ₂ /IPA/NH ₄ OH (78/20/2)	EA/MeOH/NH ₄ OH (68/28/4)

United Chemical Technologies would like to thank Dr. Leon Glass of Lab Specialties for his efforts in developing these non chlorinated mixtures.

REAGENTS

Acetic Acid, Glacial (CH_3COOH): 17.4 M

Ammonium Hydroxide (NH_4OH): concentrated (14.8 M)

β -Glucuronidase: lyophilized powder from limpets (*Patella Vulgatta*)

Dimethylformamide (DMF): silylation grade

Hydrochloric Acid (HCl): concentrated (12.1 M)

N, O-Bis(trimethylsilyl)trifluoroacetamide (BSTFA) with 1% trimethylchlorosilane (TMCS)

Pentafluoropropionic Acid Anhydride (PFAA or PFPA)

Phosphoric Acid (H_3PO_4): concentrated (14.7 M)

Sodium Acetate Trihydrate ($\text{NaCH}_3\text{COO}\cdot 3\text{H}_2\text{O}$): F.W. 136.08

Sodium Borate Decahydrate ($\text{Na}_2\text{B}_4\text{O}_7\cdot 10\text{H}_2\text{O}$): F.W. 381.37

Sodium Hydroxide, (NaOH): F.W. 40.00

Sodium Phosphate Dibasic, Anhydrous (Na_2HPO_4): F.W. 141.96

Sodium Phosphate Monobasic, Monohydrate ($\text{NaH}_2\text{PO}_4\cdot \text{H}_2\text{O}$): F.W. 137.99

NOTES:

Storage of organics in some plastic containers may lead to plasticizer contamination of the solvent or solvent mixture, this may interfering with analyte quantitation.

Good laboratory practice dictates all who handle or are potentially exposed to reagents, solvents and solutions used or stored in the laboratory should familiarize themselves with manufacturer's recommendations for chemical storage, use and handling, and should also familiarize themselves with an appropriate Material Safety Data Sheet (MSDS).

HOW TO PREPARE SOLUTIONS AND BUFFERS

Acetic Acid, 1.0 M:

To 400 mL DI H₂O add 28.6 mL glacial acetic acid. Dilute to 500 mL with DI H₂O.

Storage: 25° C in glass or plastic

Stability: 6 months

Acetic Acid, 100 mM:

Dilute 40 mL 1.0 M acetic acid to 400 mL with DI H₂O. Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Acetate Buffer, 100 mM (pH 4.5):

Dissolve 2.93 g sodium acetate trihydrate in 400 mL DI H₂O; add 1.62 mL glacial acetic acid. Dilute to 500 mL with DI H₂O. Mix. Adjust pH to 4.5 ± 0.1 with 100 mM sodium acetate or 100 mM acetic acid.

Storage: 25°C in glass or plastic

Stability: 6 months; Inspect daily with use for contamination.

Acetate Buffer, 1.0 M (pH 5.0):

Dissolve 42.9 g sodium acetate trihydrate in 400 mL DI H₂O; add 10.4 mL glacial acetic acid. Dilute to 500 mL with DI H₂O. Mix.

Adjust pH to 5.0 ± 0.1 with 1.0 M sodium acetate or 1.0 M acetic acid.

Storage: 25° C in glass or plastic

Stability: 6 months; Inspect daily with use for contamination.

Acetate Buffer, 100 mM (pH 5.0):

Dilute 40 mL 1.0 M acetate buffer to 400 mL with DI H₂O. Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Ammonium Hydroxide, 7.4 M:

To 50 mL DI H₂O add 50 mL concentrated NH₄OH. Mix.

Storage: 25° C in glass or fluoropolymer plastic

Stability: Storage condition dependent

β-Glucuronidase, *Patella Vulgata*, 5,000 Fishman units/ml:

Dissolve 100,000 Fishman units lyophilized powder with 20 mL acetate buffer, 100 mM (pH 5.0).

Storage: 5° C in plastic

Stability: Several days; prepare daily for best results.

Hydrochloric Acid, 100 mM:

To 400 mL DI H₂O add 4.2 mL concentrated HCl. Dilute to 500 mL with DI H₂O. Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Methanol/Ammonium Hydroxide (98/2):

To 98 mL CH₃OH add 2 mL concentrated NH₄OH. Mix.

Storage: 25°C in glass or fluoropolymer plastic

Stability: 1 day

HOW TO PREPARE SOLUTIONS AND BUFFERS (continued)

Methylene Chloride/Isopropanol/Ammonium Hydroxide (78/20/2):

To 40 mL IPA add 4 mL concentrated NH_4OH . Mix. Add 156 mL CH_2Cl_2 . Mix.

Storage: 25°C in glass or fluoropolymer plastic

Stability: 1 day

Phosphate Buffer, 100 mM, pH 6.0:

Dissolve 1.70 g Na_2HPO_4 and 12.14 g $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ in 800 mL DI H_2O .

Dilute to 1000 mL using DI H_2O . Mix. Adjust pH to 6.0 ± 0.1 with 100 mM monobasic sodium phosphate (lowers pH) or 100 mM dibasic sodium phosphate (raises pH).

Storage: 5° C in glass

Stability: 1 month; Inspect daily with use for contamination.

Phosphoric Acid, 500 mM:

To 400 mL DI H_2O and 17.0 mL concentrated phosphoric acid. Dilute to 500 mL with DI H_2O . Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Sodium Acetate, 1.0 M:

Dissolve 13.6 g sodium acetate trihydrate in 90 mL DI H_2O . Dilute to 100 mL with DI H_2O . Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Sodium Acetate, 100 mM:

Dilute 10 mL 1.0 M sodium acetate to 100 mL with DI H_2O . Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Sodium Borate, 100 mM:

Dissolve 3.81 g $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10 \text{H}_2\text{O}$ in 90 mL DI H_2O .

Dilute to 100 mL with DI H_2O . Mix.

Storage: 25°C in glass or plastic

Stability: 6 months

Sodium Phosphate Dibasic, 100 mM:

Dissolve 2.84 g Na_2HPO_4 in 160 mL DI H_2O . Dilute to 200 mL using DI H_2O . Mix.

Storage: 5° C in glass

Stability: 1 month; Inspect daily with use for contamination.

Sodium Phosphate, Monobasic, 100 mM:

Dissolve 2.76 g $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ in 160 mL DI H_2O . Dilute to 200 mL with DI H_2O . Mix.

Storage: 5° C in glass

Stability: 1 month; Inspect daily with use for contamination.

Sulfuric Acid, 100 mM:

To 400 mL DI H_2O add 5.6 mL concentrated H_2SO_4 . Dilute to 500 mL with DI H_2O . Mix.

Storage: 25° C in glass or plastic

Stability: 6 months

Extraction Hints

- * Verify sample application pH. Analytes that are not in their proper form (i.e., neutral or charged), will not effectively bind to the sorbent and may result in erratic recoveries.
- * Do not allow the sorbent to dry between conditioning steps or before sample application. To insure properly solvated columns, apply each solvent immediately after the previous solvent. Improperly conditioned cartridges may lead to erratic recoveries.
- * Prior to elution, fully dried cartridges will insure optimal analyte recovery. To confirm column dryness, press the sides of the cartridge at the sorbent level at full vacuum. Columns should feel ambient temperature, not cool. If the column feels cool, water is probably present. Dry the column further.
- * Always use fresh NH_4OH when preparing basic elution solvents. Proper elution pH (11-12) is critical to achieving optimal recovery of basic drugs with high pK_a 's (i.e., amphetamines, some tricyclics, morphine). NH_4OH rapidly loses its strength when exposed to air. Weak NH_4OH may lead to erratic recoveries.
- * NH_4OH is more soluble in IPA than CH_2Cl_2 . To insure complete mixing of eluate solvents, add NH_4OH and IPA, then add CH_2Cl_2 .
- * Some drugs are heat labile and will degrade if overheated. Closely monitor elution dry down to prevent loss of analyte.
- * Always condition the column with the strongest solution the column will see to insure the cleanest extraction of your eluate.
- * Solvent quantities for RSV methods are suggested and might be further reduced to meet particular laboratory needs.

AMPHETAMINES IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine add internal standard(s)* and 1 mL of 0.1M phosphate buffer (pH 6.0).
Mix/vortex.
Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 100 mM acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE AMPHETAMINES

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minutes.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. CONCENTRATE ELUATE

Add 30 µL silylation grade DMF to eluate.
Evaporate to 30 µL at < 40°C

7. FLUOROACYLATE WITH PFPA (PFAA)

Add 50 µL PFPA (PFAA). Overlay with N₂ and cap.
React 20 minutes at 70°C. Evaporate to dryness at < 40°C
Reconstitute with 100 µL ethyl acetate.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Analyte (PFA)	Primary Ion**	Secondary	Tertiary
Amphetamine	190	91	118
D ₅ Amphetamine*	194	91	123
Methamphetamine	204	118 or 91	160
D ₅ Methamphetamine*	208	119 or 92	163

*Suggested internal standards for GC/MS: D₅ Amphetamine and D₅ Methamphetamine

**Quantification Ion

**ANABOLIC STEROIDS IN URINE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1 of 2

- 1. PREPARE SAMPLE- β -GLUCURONIDASE HYDROLYSIS**
To 5 mL of urine add internal standard(s)* and 2 mL of β -glucuronidase.
 β -glucuronidase: 5,000 F units/mL Patella Vulgata in 0.1 M acetate buffer (pH 5.0).
Mix/vortex. Hydrolyze for 3 hours at 65° C. Cool before proceeding.
Centrifuge for 10 minutes at 2000 rpm and discard pellet.
Adjust sample pH to 6.0 \pm 0.5 with approximately 700 μ l of 1.0 N NaOH.
- 2. PREPARE CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
- 3. APPLY SAMPLE**
Load at 1 to 2 mL/minute.
- 4. WASH COLUMN**
1 x 3 mL 10% (v/v) CH₃OH in DI H₂O; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 1 mL hexane or hexane/ethyl acetate (50/50); aspirate.
- 5. ELUTE ANABOLIC STEROIDS (Choose a, b, c or d)**
 - a. 1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).
 - b. 1 x 3 mL CH₂Cl₂/IPA (80/20)
 - c. 1 x 3 mL ethyl acetate
 - d. 1 x 3 mL CH₃OH
- 6. DRY ELUATE**
Evaporate to dryness at < 40° C.
- 7. DERIVATIZE**
Add 50 μ l ethyl acetate and 50 μ l MSTFA (with 3% trimethylsilyliodide).
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate MSTFA solution.

**ANABOLIC STEROIDS IN URINE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

8. QUANTITATE

Inject 1 to 2 µL onto chromatograph.
Monitor the following ions (GC/MS):

Compound	Primary*	Secondary	Tertiary	OTHER
Testosterone-TMS	432	301	209	
19-Noretiocholanone-TMS	405	315	225	
Oxymethalone	640	552	462	370,143
Dehydroepiandrosterone-2TMS	432	327	297	
10-Nortestosterone-2TMS	418	287	194	
Oxymethalone Metabolite #1	640	552	462	143
Oxymethalone Metabolite #2	625	462	370	143
11-β-Hydroxyandosterone	522	417	158	
Methandienone	409	313	281	
19-Norandosterone-2TMS	405	315	225	
Alpha-Hydroxyetiocholanone	504	417		
17-α-Epitestosterone-TMS	432	341	327	209
Stanozolol	472	381	342	149

* Quantitation Ion

SOURCE- UCT Internal Publication

BARBITURATES IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine add internal standard(s)* and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.

5. ELUTE BARBITURATES

1 x 3 mL hexane/ethyl acetate (50/50); collect eluate at 1 to 2mL / minute.

6. DRY ELUATE

Evaporate to dryness at < 40° C.
Reconstitute with 100 µL ethyl acetate.

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Drug	Primary Ion**	Secondary Ion	Tertiary Ion
Amobarbital:	156	141	157
Butobarbital:	156	141	157
Butalbital:	168	167	181
Hexobarbital*	221	157	236
Pentobarbital	156	141	197
Phenobarbital	204	232	117
Secobarbital	168	167	195
Thiopental:	172	157	173

* Suggested internal standard for GC/MS: Hexobarbital or a Deuterated Barbiturate analog.

** Quantitation Ion

BARBITURATES IN URINE FOR GC OR GCMS CONFIRMATIONS USING A 80 mg CLEAN SCREEN® REDUCED SOLVENT VOLUME EXTRACTION COLUMNS

(Part # CSDAUA83 without tips or with CLEAN-THRU® Tips CCDAUA83)

1. PREPARE SAMPLE

To 1 mL of urine add internal standard(s)* and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 0.5 mL CH₃OH; aspirate.
1 x 0.5 mL DI H₂O; aspirate.
1 x 0.25 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 0.5 mL DI H₂O; aspirate.
1 x 0.5 mL 0.1 M acetic acid; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 0.1 mL hexane; aspirate.

5. ELUTE BARBITURATES

2 X 0.75 mL hexane/ethyl acetate (50/50)

6. DRY ELUATE

Evaporate to dryness at < 40° C.
Reconstitute with 100 µL ethyl acetate.

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Drug	Primary Ion**	Secondary Ion	Tertiary Ion
Amobarbital:	156	141	157
Butobarbital:	156	141	157
Butalbital:	168	167	181
Hexobarbital*	221	157	236
Pentobarbital	156	141	197
Phenobarbital	204	232	117
Secobarbital	168	167	195
Thiopental	172	157	173

* Suggested internal standard for GC/MS: Hexobarbital or a Deuterated Barbiturate analog.

** Quantitation Ion

BASIC DRUGS FOR HPLC ANALYSIS USING 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 5 mL of urine add internal standard(s) and 2 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
1 x 3 mL methanol; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE BASES

1 x 2 mL CH₃OH/NH₄OH (98/2)
Collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily.

6. EXTRACT

To eluate add 2.0 mL DI H₂O and 500 µl methylene chloride.
Mix/vortex.
Centrifuge. At 2,000 RPM for 10 minutes
Transfer organic (lower) layer to a clean test tube.

7. EVAPORATE

Evaporate to dryness at < 40° C

8. QUANTITATE

Reconstitute in mobile phase and inject onto the HPLC.

SOURCE: UCT Internal Publication

BENZODIAZEPINES IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1 of 2

1. PREPARE SAMPLE- β-GLUCURONIDASE HYDROLYSIS

To 2 mL of urine add internal standard(s)* and 1 mL of β-glucuronidase solution. β-glucuronidase solution contains 5,000 F units/mL *Patella Vulgata* in 0.1 M acetate buffer (pH=5.0). Mix/vortex. Hydrolyze for 3 hours at 65° C. Centrifuge for 10 minutes at 2000 rpm and discard pellet. Cool before proceeding.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0)
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 20% acetonitrile in 0.1 M phosphate buffer (pH 6.0); aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.

5. ELUTE BENZODIAZEPINES

1 x 3 mL ethyl acetate; collect eluate at 1 to 2 mL/minute.

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 50 μL ethyl acetate and 50 μL BSTFA (with 1% TMCS).
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.

**BENZODIAZEPINES IN URINE FOR
GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Generic Name	Trade Name	Primary Ion**	Secondary	Tertiary
Alprazolam	Xanax®	308	279	204
aHydroxyalprazolam-TMS		381	396	383
Chlordiazepoxide	Librium®	282	283	284
Clonazepam	Clonopin®	387	352	306
Diazepam	Valium®	256	283	221
Nordiazepam-TMS		341	342	343
Desalkylflurazepam-TMS		359	341	245
Hydroxyethylflurazepam-TMS		288	360	389
Lorazepam-TMS	Ativan®	429	430	347
Oxazepam-TMS	Serax®	429	430	313
Prazepam*		269	241	324
Temazepam-TMS	Restoril®	343	283	257
Triazolam	Halcion®	313	314	342
aHydroxytriazolam-TMS		415	417	430

* Suggested internal standard for GC/MS: Prazepam or D₅-Oxazepam

** Quantitation ion

Note: Flurazepam does not extract under these conditions; however metabolites such as desalkylflurazepam and hydroxyethylflurazepam will extract with high recovery.

Reference- UCT Internal Publication

**BENZODIAZEPINES IN SERUM OR PLASMA
FOR HPLC ANALYSIS:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

- 1. PREPARE SAMPLE**
To 1 mL of serum add internal standard and 1.0 mL of 0.1M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 0.1M monobasic or dibasic sodium phosphate.
- 2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1M phosphate buffer (pH 6.0)
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
- 3. APPLY SAMPLE**
Load at 1 mL/minute.
- 4. WASH COLUMN**
1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 20% acetonitrile in 0.1 M phosphate buffer (pH 6.0); aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.
- 5. ELUTE BENZODIAZEPINES**
1 x 3 mL ethyl acetate; collect eluate at 1 to 2 mL/minute.
- 6. DRY ELUATE**
Evaporate to dryness at < 40° C
- 7. RECONSTITUTE**
Reconstitute in mobile phase.
- 8. QUANTITATE**
Inject sample onto HPLC.

Reference- UCT Internal Publication

**BENZODIAZEPINES (SERUM) BY HPLC
100 mg CLEAN-UP® C2 EXTRACTION COLUMN**

(Part # CEC02111)

I. REAGENTS:

1. SATURATED SODIUM BORATE BUFFER, pH 9.5

Dissolve 10 grams of sodium borate in 950 mL of D. I. H₂O.
Adjust pH to 9.5 with 10 N sodium hydroxide. q.s. to 1000 ml.

2. POTASSIUM PHOSPHATE BUFFER, pH 6.0 (0.01M)

Dissolve 1.36 grams of monobasic potassium phosphate in 950 mL of D. I. H₂O.
Adjust pH to 6.0 with 1 N sodium hydroxide. q.s. to 1000 ml.

3. MOBILE PHASE:

18.5% acetonitrile
26.5% methanol
55.0% potassium phosphate buffer, pH 6.0 (0.01M)

II. EXTRACTION METHOD:

1. PREPARE SAMPLE:

Add 200 µL of sodium borate buffer to 0.5 mL of serum. Vortex.

2. CONDITIONING CLEAN-UP® EXTRACTION COLUMN

2 x 1 mL methanol
2 x 1 mL D.I. H₂O

DO NOT LET COLUMN DRY OUT
(vacuum settings should be 3-5 inches of Hg.)

3. APPLY SAMPLE TO COLUMN:

Add buffered sample to top of column.
Pull through at a flow of 1-2 ml/minute.

4. WASH COLUMN:

1 x 1 mL D.I.H₂O
Dry column 2-3 minutes under vacuum (15-20 inches of Hg).

5. ELUTE BENZODIAZEPINES:

2 x 0.5 mL methanol

6. DRY ELUATE:

Dry eluate under nitrogen and minimal heat.
Reconstitute with 100 µL of mobile phase.

7. QUANTITATE

Inject 15-20 µL onto HPLC.

Reference - UCT Internal Publication

BETA BLOCKERS IN URINE FOR GC/MS CONFIRMATIONS USING 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **CONDITION CLEAN SCREEN® EXTRACTION COLUMN**

1 x 3 mL methanol; aspirate.

1 x 3 mL deionized water; aspirate.

1 x 3 mL 0.1 M Acetate Buffer (pH 4.7)

NOTE: Aspirate at ≤ 3 inches Hg to prevent sorbent drying.

2. **Apply sample**

Take one milliliter of urine and add 2 mLs of 0.1M Acetate Buffer at pH 4.7

Load at 1 to 2 mL/ minute.

3. **Wash column**

2 X 1 mL Acetone/ Methanol (1:1) aspirate.

Dry column (5 minutes at ≥ 10 inches Hg).

4. **Elute Beta Blockers**

1 x 1 mL Dichloromethane/ Isopropanol and Ammonium Hydroxide (82/16/2)

Collect the eluate by gravity

NOTE: Prepare elution solvent fresh daily.

Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

5. **Dry Eluate**

Evaporate to dryness at $\leq 40^\circ$ C.

6. **Derivatize**

Derivatization Solution- Methaneboronic acid at 5mg/mL was prepared in dry ethyl acetate (use molecular sieve)

Store this solution at -20° C (freezer conditions) until use.

Reaction Mixture

Add 100 μ L of the Methaneboronic acid solution (see above)

Mix/vortex.

React 15 minutes at 70° C. Remove from heat source to cool.

NOTE: Do not evaporate this solution.

7. **Analysis**

Inject 1 to 2 μ L sample

Reference:

Branum G, Sweeney S, Palmeri A, Haines L and Huber C

The Feasibility of the Detection and Quantitation of β Adrenergic Blockers By Solid Phase Extraction and Subsequent Derivatization with Methaneboronic Acid.

Journal of Analytical Toxicology 22: 135-141 (1998)

BETA AGONISTS IN URINE FOR GC/MS CONFIRMATIONS USING 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

- 1. CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL methanol; aspirate.
1 x 3 mL deionized water; aspirate.
1 x 3 mL 0.1 M Acetate Buffer (pH 4.7)
NOTE: Aspirate at ≤ 3 inches Hg to prevent sorbent drying.
- 2. Apply sample**
Take one milliliter of urine and add 2 mLs of 0.1M Acetate Buffer at pH 4.7
Load at 1 to 2 mL/ minute.
- 3. Wash column**
2 X 1 mL Acetone/ Methanol (1:1) aspirate.
Dry column (5 minutes at ≥ 10 inches Hg).
- 4. Elute Beta Blockers**
1 x 1 mL Dichloromethane/ Isopropanol and Ammonium Hydroxide (82/16/2)
Collect the eluate by gravity feed
NOTE: Prepare elution solvent fresh daily.
Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).
- 5. Dry Eluate**
Evaporate to dryness at $\leq 40^\circ$ C.
- 6. Derivatize**
Derivatization Solution-
Methaneboronic acid at 5mg/mL was prepared in dry ethyl acetate (use molecular sieve)
Store this solution at -20° C (freezer conditions) until use.
Add 100 μ L of the Methaneboronic acid solution (see above)
Mix/vortex.
React 15 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate this solution.
- 7. Analysis**
Inject 1 to 2 μ L sample (derivatized solution)

CARBOXY THC IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSTHC020 without Tips or ZCTHC020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE - BASE HYDROLYSIS OF GLUCURONIDES

To 2 mL of urine add internal standard* and 100 µL of 10 N NaOH.
Mix/vortex.
Hydrolyze for 20 minutes at 60° C. Cool before proceeding.
Adjust sample pH to 3.5 ± 0.5 with 1.0 ml of glacial acetic acid.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL H₂O; aspirate.
1 x 1 mL 0.1M HCl; aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 ml/minute.

4. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 0.1M HCl/acetonitrile (70/30); aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 200 µL hexane; aspirate.

5. ELUTE CARBOXY THC

1 x 3 mL hexane/ethyl acetate (50/50) Collect eluate at 1 to 2 ml / minute

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 50 µL ethyl acetate and 50 µL BSTFA (with 1% TMCS).
Overlay with Nitrogen and cap.
Mix/vortex.
React 20 minutes at 70° C.
Remove from heat source to cool.
NOTE: Do not evaporate BSTFA.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

ANALYTE (TMS)	Primary Ion**	Secondary	Tertiary
Carboxy-Δ ⁹ THC	371	473	488
D ₃ Carboxy-Δ ⁹ THC*	374	476	491

*Suggested internal standard for GC/MS: D₃Carboxy-Δ⁹THC

**Quantitation ion.

Reference- UCT Internal Publication

**CARBOXY THC IN URINE MANUAL METHOD
FOR GC OR GC/MS CONFIRMATIONS USING:
80 mg CLEAN SCREEN®
REDUCED SOLVENT VOLUME EXTRACTION COLUMN**

(Part # CSDAUA83 without Tips or CCDAUA83 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE - BASE HYDROLYSIS OF GLUCURONIDES

To 1 mL of urine add internal standard* and 100 µL of 10 N NaOH.
Mix/vortex.
Hydrolyze for 20 minutes at 60° C. Cool before proceeding.
Adjust sample pH to 3.5 ± 0.5 with about 400 µL of glacial acetic acid.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 500 µL Hexane/ Ethyl Acetate (1:1)
1 x 500 µL CH₃OH; aspirate.
1 x 500 µL H₂O; aspirate.
1 x 250 µL 0.1M HCl; aspirate.

NOTE: Use gravity flow or minimal vacuum to condition the column

3. APPLY SAMPLE

Load at 1 to 2 ml/minute.

4. WASH COLUMN

1 x 500 µL DI H₂O; aspirate.
1 x 500 µL 0.1M HCl/acetonitrile (70/30); aspirate.
Dry column (3 minutes at > 10 inches Hg).
1 x 100 µL hexane; aspirate.

5. ELUTE CARBOXY THC

2 x 750 µL hexane/ethyl acetate (50/50).

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 50 µL ethyl acetate and 50 µL BSTFA (with 1% TMCS).
Overlay with Nitrogen and cap.
Mix/vortex.
React 20 minutes at 70° C.
Remove from heat source to cool.

NOTE: Do not evaporate BSTFA.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

ANALYTE (TMS)	Primary Ion**	Secondary	Tertiary
Carboxy-Δ ⁹ THC	371	473	488
D ₃ Carboxy-Δ ⁹ THC*	374	476	491

*Suggested internal standard for GC/MS: D₃Carboxy-Δ⁹ THC

**Quantitation ion.

CARISOPRODOL AND MEPROBAMATE FOR GC OR GC/MS CONFORMATIONS USING A 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **PREPARE SAMPLE**
To 0.5 mLs of urine sample add internal standard and 2mLs of 0.1M pH 6.0 Phosphate buffer
2. **CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL methanol; aspirate.
1 x 3 mL deionized water; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0)
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
3. **APPLY SAMPLE**
Load at 1 to 2 mL/ minute.
4. **WASH COLUMN**
1 x 4 mL deionized water; aspirate.
1 x 2 mL 0.1M acetic acid; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 3 mL hexane; aspirate.
5. **ELUTE CARISOPRODOL/MEPROBAMATE**
1 x 3 mL hexane/ethyl acetate (50:50)
Collect eluate at 1 to 2 mL/minute.
6. **DRY ELUATE**
Evaporate to dryness at < 40° C.
7. **QUANTITATE**
Reconstitute with 100 µL ethyl acetate
Inject 1 to 2 µL sample on gas chromatograph
For MSD monitor the following ions:

Compound	Primary Ion**	Secondary	Tertiary
Carisoprodol	158	104	245
Meprobamate	83	114	144
Hexobarbital *	221	157	81

*Suggested internal standard for GC/MS: Hexobarbital

**Quantitation ion

CLONAZEPAM & 7-AMINOCLONAZEPAM IN URINE FOR GC/MS CONFIRMATIONS USING 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **Prepare sample - β -Glucuronidase hydrolysis**
To 2 mL of urine add internal standard(s)* and 1 mL of β -glucuronidase solution.
 β -glucuronidase solution contains 5,000 F units/mL *Patella Vulgata* in 0.1 M acetate buffer (pH 5.0).
Mix/vortex. Hydrolyze for 3 hours at 65° C
Cool before proceeding.
2. **CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL methanol; aspirate.
1 x 3 mL deionized water; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0)
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
3. **Apply sample**
Load at 1 to 2 mL/ minute.
4. **Wash column**
1 x 2 mL deionized water; aspirate.
1 x 2 mL 20% acetonitrile in 0.1 M phosphate buffer (pH 6.0); aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.
5. **Elute Clonazepam / 7-Aminoclonazepam**
1 x 3 mL ethyl acetate with 2% NH₄OH: collect eluate at 1 to 2 mL/minute.
NOTE: Prepare fresh daily.
6. **Dry Eluate**
Evaporate to dryness at < 40° C.
7. **Derivatize**
Add 50 μ L ethyl acetate and 50 μ L MTBSTFA (with 1% TBDMCS).
Mix/vortex.
React 20 minutes at 90° C. Remove from heat source to cool.
NOTE: Do not evaporate MTBSTFA solution.
8. **Analysis**
Inject 1 to 2 μ L sample
For MSD monitor the following ions:

Compound (TBDMS)	Primary Ion**	Secondary	Tertiary
Clonazepam	372	374	326
7-Aminoclonazepam	342	344	399
Oxazepam-D ₅ *	462	464	463

*Suggested internal standard for GC/MS: Oxazepam-D₅

**Quantitation ion

COCAINE AND BENZOYLECGONINE IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **PREPARE SAMPLE**
To 2 mL urine add internal standard(s)* and 1 mL of 0.1 M phosphate buffer (pH 6.0). Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.
2. **CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
3. **APPLY SAMPLE**
Load at 1 to 2 mL/minute.
4. **WASH COLUMN**
1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 0.1 M HCl; aspirate.
1 x 3 mL methanol; aspirate.
Dry column (5 minutes at > 10 inches Hg).
5. **ELUTE COCAINE AND BENZOYLECGONINE**
1 x 3 mL Methylene Chloride/Isopropanol/Ammonium Hydroxide (78/20/2)
Collect eluate at 1 to 2 mL/minutes.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).
6. **DRY ELUATE**
Evaporate to dryness at < 40° C
7. **DERIVATIZE**
Add 50 µL ethyl acetate and 50 µL BSTFA (with 1% TMCS).
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.
8. **QUANTITATE**
Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Cocaine	182	198	303
D ₃ -Cocaine*	185	201	306
Benzoylcegonine-TMS	240	256	361
D ₃ -Benzoylcegonine-TMS*	243	259	364

* Suggested internal standards for GC/MS: D₃-Cocaine, D₃-Benzoylcegonine

**Quantitation Ion

**COCAINE AND BENZOYLECGONINE
IN SERUM, PLASMA, OR WHOLE BLOOD
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 1 mL of sample (Serum, Plasma or Whole Blood) add internal standard and 4 mL of DI H₂O.
Mix/vortex and let stand 5 minutes.
Centrifuge for 10 minutes at 2,000 rpm and discard pellet.
Add 2 mL of 0.1 M phosphate buffer (pH 6.0). Mix/vortex.
Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 0.1 M HCl; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE COCAINE AND BENZOYLECGONINE

1 x 3 mL CH₂Cl₂ /IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 50 µL ethyl acetate and 50 mL BSTFA (with 1% TMCS).
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Cocaine	182	198	303
D ₃ -Cocaine*	185	201	306
Benzoyllecgonine -TMS	240	256	361
D ₃ -Benzoyllecgonine -TMS*	243	259	364

*Suggested internal standards for GC/MS: D₃-Cocaine, D₃-Benzoyllecgonine

** Quantitation Ion

COCAINE AND BENZOYLECGONINE IN SERUM, PLASMA, OR WHOLE BLOOD FOR HPLC USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 1 mL of sample (Serum, Plasma or Whole Blood) add internal standard(s) and 4 mL of DI H₂O.

Whole Blood: Mix/vortex and let stand 5 minutes.

Centrifuge for 10 minutes at 2000 rpm and discard pellet.

Add 2 mL of 0.1M phosphate buffer (pH 6.0). Mix/vortex.

Sample pH should be 6.0 ± 0.5.

Adjust pH accordingly with 0.1M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.

1 x 3 mL DI H₂O; aspirate.

1 x 1 mL 0.1M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.

1 x 2 mL 0.1 M HCl; aspirate.

1 x 3 mL CH₃OH; aspirate.

Dry column (5 minutes at > 10 inches Hg).

5A*. ELUTE COCAINE AND BENZOYLECGONINE

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

5B*. ELUTE COCAINE AND BENZOYLECGONINE

1 x 2 mL CH₃OH/NH₄OH (98/2); collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent daily.

Add 3 mL DI H₂O and 500 µL CH₂Cl₂ to eluate.

Mix / vortex 10 seconds. Centrifuge if necessary to separate layers.

Aspirate and discard aqueous (upper) layer.

6. CONCENTRATE

Evaporate to dryness at < 40° C.

Reconstitute in mobile phase for injection into HPLC.

* Choose either 5A or 5B

COCAINE AND ITS METABOLITES FROM MECONIUM FOR GC OR GC/MS ANALYSIS USING CLEAN SCREEN® EXTRACTION COLUMNS

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

Vortex 0.5 -1 g meconium and 2 mL of CH₃OH.
Centrifuge and transfer the supernatant to a clean tube.
To each tube add 3 mL 0.1M phosphate buffer (pH 6.0), internal standard and vortex.
Matrix must be more aqueous than organic for good retention to occur.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

2 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 3 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute. Allow to dry

4. WASH COLUMN

1 x 1 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M HCl; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE COCAINE AND METABOLITES

1 x 3 mL CH₂Cl₂/IPA/NH₂OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. EVAPORATE

Evaporate the elution solvent to dryness without heating.

7. DERIVATIZE

Add 50 µL ethyl acetate and 50 µL BSTFA (with 1% TMCS).
Overlay with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Cocaine	182	198	303
D ₃ -Cocaine*	185	201	306
Benzoylcegonine-TMS	240	256	361
D ₃ -Benzoylcegonine-TMS*	243	259	364

* Suggested internal standards for GC/MS

** Quantitation ion

DHEA, TESTOSTERONE, AND EPITESTOSTERONE IN URINE FOR GC OR GC/MS ANALYSIS USING: 200 mg CLEAN THRU® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **PREPARE SAMPLE**
Pipette five milliliters of urine into 16 x 100 mm borosilicate glass test tubes.
Add internal standard*, adjust sample pH to 5.5- 6.5 using concentrated sodium phosphate monobasic or dibasic. Mix sample.
Centrifuge samples at 3000 rpm for 5 min.
2. **CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 3 mL 0.1M pH 6.0 phosphate.
3. **APPLY SAMPLE**
Pour supernatant onto column.
Allow to flow via gravity.
4. **WASH COLUMN**
1 x 3 mL DI H₂O; aspirate.
Dry column (10 minutes at > 10 mm Hg).
5. **ELUTE STEROIDS**
1 x 3 mL of methanol.
6. **ENZYMATIC HYDROLYSIS**
Eluates were dried under a stream of nitrogen after which 2 mL of 0.2M phosphate buffer (pH 7.0) and 250 units of β-glucuronidase were added, vortex mixed, and allowed to incubate at 50° C for 1 hour. Samples were then cooled, and the pH was adjusted to 10-11 using a 1:1 mixture of NaHCO₃/Na₂ CO₃.
7. **ADDITION CLEAN-UP®**
Add 5 mL of n-butyl chloride to each sample. The tubes were capped and shaken vigorously for 10 minutes and then centrifuged at 3000 rpm for 5 min. the organic layer was transferred to clean test tubes and dried under a stream of nitrogen. Place dried sample in a dessicator and further dry under vacuum for 30 minutes.
8. **DERIVATIZE**
Add 50 μL of MSTFA/NH₄I/dithioerythritol (1000:2:5, V/W/W) and incubate at 70° C for 20 min.
Centrifuge sample at 3000 rpm for 1 min. and transfer directly to gc injector vials.
9. **QUANTITATE**
Inject 1 to 2 μL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary Ion **	Secondary Ion
Testosterone	432	417
Epitestosterone	432	417
DHEA	432	417
16 a Hydroxytestosterone*	520	

* Suggested internal standard at (20 ng/mL)

** Quantitation ion

SOURCE- UCT Internal Publication

FENTANYL AND ANALOGUES IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 5 mL of sample add internal standard* and 2 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE FENTANYLS

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. CONCENTRATE

Evaporate to dryness at < 40° C.
Reconstitute with 50 µl ethyl acetate.

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Fentanyl	245	146	189
D ₅ -Fentanyl*	250	151	194
Alpha-Methylfentanyl	259	203	146
Para-Fluorofentanyl	263	164	207
3-Methylfentanyl	259	160	203
Thienfentanyl	245	146	189
Sufentanil	289	140	
Carfentanil	303	187	
Lofentanil	317	201	289
Alfentanil	289	268	222

* Suggested internal standard for GC/MS: D₅-Fentanyl

** Quantitation ion

SOURCE- UCT Internal Publication working with the Philadelphia Medical Examiner's Office

FLUNITRAZEPAM AND METABOLITES IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1 of 2

1. PREPARE SAMPLE- β -GLUCURONIDASE HYDROLYSIS

To 2 mL of urine add internal standard(s)* and 1 mL of β -glucuronidase solution.
 β -glucuronidase solution contains 5,000 F units/mL *Patella Vulgata* in 0.1 M acetate buffer (pH=5.0). Mix/vortex. Hydrolyze for 3 hours at 65° C.
Centrifuge for 10 minutes at 2000 rpm and discard pellet.
Cool before proceeding.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0)
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 20% acetonitrile in 0.1 M phosphate buffer (pH 6.0); aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.

5. ELUTE FLUNITRAZEPAM, 7-AMINOFLUNITRAZEPAM AND DESMETHYLFLUNITRAZEPAM

1 x 3 mL ethyl acetate with 2% NH₄OH; collect eluate at 1 to 2 mL/minute.
Prepare fresh daily.

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 50 μ L ethyl acetate and 50 μ L MTBSTFA (with 1% TBDMCS).
Overlay with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.

NOTE: Do not evaporate MTBSTFA solution.

**FLUNITRAZEPAM AND METABOLITES IN URINE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

2 of 2

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.

For MSD monitor the following ions:

Compound	Primary Ion**	Secondary	Tertiary
Flunitrazepam	312	286	266
7-Aminoflunitrazepam	283	255	254
Desmethyflunitrazepam	356	357	310
Oxazepam-D ₅ ISTD*	462	464	463

* Suggested internal standard for GC/MS: Oxazepam-D₅ ISTD

** Quantitation ion

GC CONDITIONS

Column-DB5 or equivalent capillary column (15 meters x 0.25mm ID x 0.25µ Film)

Injector Temperature = 250°C

INJECTOR Splitless mode

Temperature Program 180°C to 275°C at 10°C/min then

275°C to 300°C at 25°C/min

Reference- UCT Internal Publication

GABAPENTIN IN SERUM, PLASMA, OR WHOLE BLOOD FOR GC OR GC/MS ANALYSIS USING: 100 mg CLEAN UP® C18 EXTRACTION COLUMNS

(Part # CUC18111)

1. PREPARE SAMPLE

500 µL of sample, calibrator, or control was placed into a 10 x 25 mm disposable glass test tube and 25 µL of internal standard* (5.0mg/l) was added. Vortex tube. Add 500 µL of 20% acetic acid and vortex tube again.

2. CONDITION CLEAN UP® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate
1 x 3 mL DIH₂O; aspirate
1 x 3 mL 0.1N HCL aspirate

3. APPLY SAMPLE

Load at 1 to mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate
1 x 3 mL ethyl acetate
1 x 3 mL hexane
Dry column (5 minutes at > 10 inches Hg) or until column is dry

5. ELUTION

1 x 1 mL 2% NH₄OH in MeOH

6. DRY ELUATE

Evaporate to dryness at < 40°C

7. DERIVATIZATION

Add 50 µL of MTBSTFA + 1 % t-BDMCS
Cap and heat at 70° C for 30 minutes.
Remove and allow to cool

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.

INTERNAL STANDARD:

* 1-aminomethyl-1-cycloheptyl acetic acid

Reference: Carl E. Wolf II, Joseph Sady, and Alphonse Pokalis
Determination of Gabapentin in Serum using Solid Phase
Extraction and Gas-Liquid Chromatography.
Journal of Analytical Toxicology 20:498-501 (October 1996)

**A Solid Phase method for
Gamma-Hydroxybutyrate (GHB)
In Urine without conversion to Gamma-Butyrolactone (GBL)**

(Part # ZSGHB020)

Developed by:
United Chemical Technologies, Inc.

2731 Bartram Road
Bristol, PA 19007
215-781-9255
215-785-1226 fax

- 1. PREPARE SAMPLE**
To 200 μ L of urine add internal standard (GHB-D6) and 100 μ L of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex.
- 2. CONDITION CLEAN SCREEN[®] GHB EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 0.5 mL 0.1 M phosphate buffer (pH6.0); aspirate.
Note: Aspirate at < 3 inches of Hg to prevent sorbent drying.
- 3. APPLY SAMPLE**
Place test tubes into vacuum manifold for collection.
The sample loading and wash are both collected.
Decant sample onto column. Aspirate at ~1 inch Hg.
- 4. WASH COLUMN**
Add 1 mL of CH₃OH /NH₄OH (99/1) to original sample test tube; vortex.
Decant wash onto column.
Note: Aspirate at ~1 inch of Hg.
- 5. CONCENTRATE**
Remove test tubes from vacuum manifold.
Evaporate to dryness at 60° C using a stream of air or N₂.
- 6. SAMPLE CLEAN UP**
Add 200 μ L of dimethylformamide.
Add 1 mL of hexane saturated with dimethylformamide.
Mix by inversion for 5 minutes.
Centrifuge at 3000 rpm for 5 minutes.
Transfer lower dimethylformamide layer to a clean test tube.
- 7. CONCENTRATE**
Evaporate to dryness at < 50° C using a stream of air or N₂.
- 8. DERIVATIZE**
Add 100 μ L ethyl acetate and 100 μ L BSTFA (with 1% TMCS).
Mix/vortex.
- 9. QUANTITATE**
Inject 1 to 2 μ L onto gas chromatograph.
For MSD monitor the following ions:

GHB-diTMS	233**, 234, 235
GHB-D6-diTMS*	239**, 240, 241

* Suggested internal standard for GC/MS: GHB-D6-diTMS

** Quantitation ion

Blood GHB Extraction Procedure
Using:
(United Chemical Technologies, Inc. Worldwide Monitoring[®] part number ZSGHB020)

by: Mr. Jim Oeldrich
Wisconsin State Crime Lab Milwaukee, WI

- 1. PREPARE SAMPLE**
To 1 mL blood sample add internal standard and 0.5 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex.
Rock for 10 minutes
Centrifuge for 10 minutes at 2700 rpm.
- 2. CONDITION CLEAN SCREEN[®] GHB EXTRACTION COLUMN**
1 x 3 mL Methanol; aspirate.
1 x 3 mL DI Water; aspirate.
1 x 0.5 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
Note: Aspirate at less than 3 inches of Hg to prevent sorbent drying.
- 3. APPLY SAMPLE**
Place centrifuge tubes into vacuum manifold for collection.
The sample loading is collected.
Decant sample onto column. Aspirate at about 1 inch Hg.
After the sample is off the columns apply full vacuum for about 15 seconds to remove any residual blood.
- 4. ELUTE GHB**
Remove centrifuge tubes, set aside.
Place clean centrifuge tubes into vacuum manifold for collection.

1 x 2 mL of CH₃OH /NH₄OH (99/1). Aspirate at about 1 inch of Hg.
- 5. CONCENTRATE**
Remove test tubes from vacuum manifold.
Vortex the sample prior to concentrating.
Evaporate to dryness at 60° C using a stream of nitrogen.
- 6. SAMPLE CLEAN UP**
Add 200 µL of dimethylformamide.
Add 1 mL of hexane saturated with dimethylformamide.
Rock for 5 minutes.
Centrifuge at 5 minutes at 2700 rpm.
Transfer lower dimethylformamide layer to a clean test tube.
(If necessary transfer all liquid to a clean tube and allow to separate, then proceed to extract the lower layer)
- 7. CONCENTRATE**
Evaporate to dryness at 50° C using a stream of air or Nitrogen.
- 8. DERIVATIZE**
Add 25 µL ethyl acetate and 25 µL BSTFA w/1% TMCS
Mix/vortex.
Heat at 70° C for 30 minutes.
- 9. QUANTITATE**
Inject a 1 to 2 µL of the sample onto GC/MS.

To Order Call: 800-541-0559 • Fax: 215-785-1226 • Shop Online: www.unitedchem.com

**A Solid Phase method for
Gamma-Hydroxybutyrate (GHB)
In Blood, Urine, Vitreous or Tissue
without conversion to Gamma-Butyrolactone (GBL)**

(Part # ZSGHB020)

1 of 2

Developed by:
Saint Louis University Health Sciences Center
Saint Louis University Medical School
Mr. Joseph A. Crifasi, M.A., M.T., (ASCP)
Certified Toxicology Specialist, ABFT

Forensic Toxicology
6030 Helen Ave.
St. Louis, MO 63134
314-522-6410 ext. 6517
314-522-0955 fax

GHB working standard; 200 µg/mL in H₂O; prepared from Radian stock 1 mg/mL.

GHB-D₆ working internal standard; 100 µg/mL; use as supplied Radian stock (0.1 mg/mL).

Working Standard	Whole Blood	Concentration
10 µL	200 µL	10 µg/mL
25µL	200 µL	25 µg/mL
50 µL	200 µL	50µg/mL
100 µL	200 µL	100 µg/mL

- 1) Make calibration standards and pipet 200 µL of QC and unknown bloods* into appropriately labeled 1.5 mL plastic centrifuge tubes.

***ALL SAMPLES INCLUDING URINE, VITREOUS OR HOMOGENIZED TISSUES (1:4)**

- 2) Add 25 µL of internal standard.
- 3) Add 1 mL of acetone; vortex 15 seconds
- 4) Centrifuge; transfer acetone layer to culture tubes.
- 5) Evaporate extracts @ 80°C w/Nitrogen
- 6) Reconstitute the dried extracts with 200 µL of 0.1 Phosphate Buffer pH 6.0 buffer; vortex 15 seconds.
- 7) **Prepare CLEAN SCREEN® GHB** (200 mg in a 10 mL tube)
Extraction columns as follows:

- 1 x 3 mL of CH₃OH; aspirate.
- 1 x 3 mL of DIH₂O; aspirate.
- 1 x 3 mL of 0.1 M Phosphate Buffer (pH 6.0) aspirate.

NOTE: Aspirate at 3 inches of Hg or less to prevent sorbent drying.

**A Solid Phase method for
Gamma-Hydroxybutyrate (GHB)
In Blood, Urine, Vitreous or Tissue
without conversion to Gamma-Butyrolactone (GBL)**

2 of 2

8) APPLY SAMPLE

Add sample with Eppendorf pipette
Aspirate at ~1 inch Hg.

9) ELUTE GHB

Place clean test tubes into vacuum manifold
Add 1 mL of CH₃OH/NH₄OH (99/1) to original sample test tube; vortex
Decant onto column and collect extract.
Aspirate ~1 inch Hg.

10) Concentrate

Remove test tube from Vacuum Manifold.
Evaporate to dryness at 70°C using a stream of Nitrogen or air.

11) Derivatize

Add 100 µL of ethyl acetate and 100 µL of BSTFA with 1% TCMS. Mix/Vortex.
No heating is required.

12) QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

GHB-diTMS	233**, 234, 235
GHB-D6-diTMS*	239**, 240, 241

* Suggested internal standard for GC/MS: GHB-D6-diTMS

** Quantitation ion

Quality Control NOTE:

Quality control samples were prepared using drug free blood and 1mg/mL in house stock standard prepared using GHB stock from Sigma (#H-3635). A negative, low and high QC sample was prepared and stored frozen in 0.5-mL aliquots until use.

**GLYCOPYRROLATE (ROBINUL) FROM EQUINE
URINE BY LC/MS/MS USING
500 mg CLEAN UP® CCX2 EXTRACTION COLUMN**

(Part # CUCCX25Z)

1. SAMPLE PREPARATION

Buffer 5 mL of urine to pH 7.0 by adding 3 mL of 0.1M phosphate buffer (pH 7.0)
Add (12.5 ng) of mepenzolate (internal standard)
Add 5 mL of water to the sample
Vortex or shake thoroughly
Centrifuge for 5 min at 800 rpm

2. CONDITION CLEAN UP® EXTRACTION COLUMN

2 x 2.5 mL MeOH
2 x 2.5 mL DiH₂O
2 x 2.5 mL phosphate buffer (0.1M, pH 7.0)

3. APPLY SAMPLE

Decant supernatant onto SPE column
Load at 1 to 2 mL / min.

4. WASH COLUMN

5 mL of MeOH
5 mL of H₂O
Dry column (5 mins > 10 inches Hg)

5. ELUTE GLYCOPYRROLATE

1 x 4 mL MeOH/0.5M NH₄OAC buffer, pH 3.0 (95:5)

6. DRY ELUTE

Evaporate to dryness at 60°C
Reconstitute with 100 µL MeOH

7. QUANTITATE

Inject 10 µL onto HPLC

KETAMINE IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine add internal standard and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. PREPARE CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1M Acetic Acid; aspirate
1 x 3 mL Methanol aspirate
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE KETAMINE

1 x 3 mL Dichloromethane/ Isopropanol/ Ammonium Hydroxide (78/20/2)
Collect eluants at 1-2 mLs/min using minimal vacuum
NOTE: Make the elution solvent fresh daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C.
Reconstitute with 100 µL methanol

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Ketamine	180	209	152
D4-Ketamine*	184	213	156

* Suggested internal standard for GC/MS: D4-Ketamine

** Quantitation ion

SOURCE- UCT Internal Publication

**LYSERGIC ACID DIETHYLAMIDE (LSD)
IN SERUM, PLASMA, OR WHOLE BLOOD
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 1 mL of sample (serum, plasma, or whole blood) add 4 mL deionized water and internal standard.*
Mix/vortex and let stand 5 minutes.
Centrifuge for 10 minutes at 2000 rpm and discard pellet.
Add 2 mL 100 mM phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 100 mM phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg. to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 100 mM acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate. Dry column (5 minutes at > 10 inch

5. ELUTE LSD

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 20 µL ethyl acetate and 20 µL BSTFA (with 1% TMCS).
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

LSD-TMS	395**, 293, 268
D ₃ -LSD-TMS*	298**, 296, 271

* Suggested internal standard for GC/MS: D₃-LSD

** Quantitation ion

**LYSERGIC ACID DIETHYLAMIDE (LSD)
IN URINE FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **PREPARE SAMPLE**
To 5 mL of Urine add internal standard and 2 mL 100 mM phosphate buffer (pH 6.0). Mix/vortex. Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.
2. **CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 100 mM phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg. to prevent sorbent drying.
3. **APPLY SAMPLE**
Load at 1 mL/minute.
4. **WASH COLUMN**
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 100 mM acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).
5. **ELUTE LSD**
1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).
6. **DRY ELUATE**
Evaporate to dryness at < 40° C.
7. **DERIVATIZE**
Add 20 µL ethyl acetate and 20 µL BSTFA (with 1% TMCS).
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.
8. **QUANTITATE**
Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

LSD-TMS	395**, 293, 268
D ₃ -LSD-TMS*	298**, 296, 271

* Suggested internal standard for GC/MS: D₃-LSD

** Quantitation ion

MANUAL METHOD FOR IMMUNOASSAY: PRELIMINARY SCREENING IN WHOLE BLOOD USING A CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 1 mL of blood add 4 mL of H₂O (5 < pH < 7).
Mix/vortex. Let stand for 5 minutes to lyse red blood cells.
Centrifuge for 10 minutes at 2000 rpm and discard pellet.
Add 2 mL of 0.1 M phosphate buffer (pH 6.0). Mix/vortex.
Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL methanol; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O
1 x 1 mL 0.1 M acetic acid; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.

5. ELUTE ACIDIC AND NEUTRAL DRUGS

1 x 3 mL hexane/ethyl acetate (50/50)
Collect eluate at < 5 mL/minute.
Remove collection tubes

6. WASH COLUMN

1 x 3 mL methanol; aspirate.
Dry column (5 minutes at > 10 inches Hg).

7. ELUTE BASIC DRUGS

Replace collection tubes from step 5
1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.

NOTE: Elute into tubes containing the acidic and neutral drugs.
Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

8. DRY ELUATE-COMBINE ELUATES (STEP 5 & 7)

Evaporate to a volume 100 µL at < 40° C

9. RECONSTITUTE

Add 900 µL of DI H₂O (Sample volume is now its original 1.0 mL).

10. ANALYZE BY EMIT

Process according to urine drug screening protocols provided by immunassay manufacturer.

SOURCE: UCT Internal Publication

To Order Call: 800-541-0559 • Fax: 215-785-1226 • Shop Online: www.unitedchem.com

METHAQUALONE IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. **PREPARE SAMPLE**
To 2 mL of urine add internal standard and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.
2. **PREPARE CLEAN SCREEN® EXTRACTION COLUMN**
1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
3. **APPLY SAMPLE**
Load at 1 mL/minute.
4. **WASH COLUMN**
1 x 3 mL DI H₂O; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.
5. **ELUTE METHAQUALONE**
1 x 3 mL hexane/ethyl acetate (50/50); collect eluate at < 5 mL/minute.
6. **DRY ELUATE**
Evaporate to dryness at < 40° C.
Reconstitute with 100 µl ethyl acetate.
7. **QUANTITATE**
Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Methaqualone	235	250	233
Hexobarbital*	221	157	156

* Suggested internal standard for GC/MS: Hexobarbital

** Quantitation ion

SOURCE- UCT Internal Publication

METHADONE IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine add internal standard(s) and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 0.1 M Monobasic or 0.1M Dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 2 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE METHADONE

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. CONCENTRATE

Evaporate to dryness at < 40° C.
Reconstitute with 100 µL ethyl acetate.

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary
Methadone	72	223	294
Methadone D ₉ *	78	226	303

* Suggested internal standard for GC/MS: Methadone D₉

** Quantitation ion

SOURCE- UCT Internal Publication

METHYLMALONIC ACID FROM SERUM OR PLASMA FOR GC/MS ANALYSIS USING 500 mg CLEAN UP® QAX EXTRACTION COLUMN

(Part # CUQAX15Z)

- 1. PREPARE SAMPLE**
Add 100 µL of internal standard d3 MMA and 1 mL of acetonitrile to 1 mL of plasma or serum
Vortex for 20 sec.
Centrifuge for 5 min at 2000 rpm
- 2. CONDITION CLEAN UP® EXTRACTION COLUMN**
Wash with 1 x 3 mL MeOH
Wash with 1 x 5 mL DI H₂O
- 3. APPLY SAMPLE**
Decant supernatant onto SPE column
- 4. WASH COLUMN**
1 x 10 mL of deionized water
Dry with vacuum for 3 min.
1 x 5 mL of methanol
Dry with vacuum for 3 min
1 x 2 mL of MTBE
Dry with vacuum for 3 min.
- 5. ELUTE METHYLMALONIC ACID**
1 x 5 mL of 3% formic acid in MTBE, collect at 1 to 2 mL/min.
- 6. DRY ELUATE**
Dry under a stream of nitrogen at < 35°C
- 7. DERIVATIZE**
Reconstitute with 25 µL of MSTFA + 1% TMCS and 25 µL acetonitrile
Heat for 20 mins at 60°C
- 8. QUANTITATE**
Inject 1 to 2 µL onto gas chromatograph.

Compliments of Mark M. Kusmin and Gabor Kormaromy-Hiller
ARUP LABORATORIES

NICOTINE AND COTININE IN URINE OR SERUM FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine or serum add internal standard(s)* and 2 mL of 0.1M phosphate buffer (pH 6.0). Mix/vortex.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 2 mL 0.1M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL Hexane; aspirate

5. ELUTE COTININE

1 x 3 mL Hexane/Ethyl acetate (50:50), collect at 1-2 mL/min.

6. WASH COLUMN

Remove rack of collection tubes to rewash columns
1 x 3 mL MeOH, aspirate.
Dry column, 5 minutes at > 10 inches Hq)

7. ELUTE NICOTINE

Replace rack of collection tubes
1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 mL/minute.

NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

8. CONCENTRATE

Evaporate to dryness at < 40° C.
Take care not to over-heat or over-evaporate
Reconstitute with 100 µL methanol

9. QUANTITATE

Inject 1 to 2 µL onto chromatograph.
Monitor the following ions (GC/MS):

Compound	Primary*	Secondary	Tertiary
Nicotine	84	133	162
Cotinine	98	119	176

* Quantitation Ion

** D₃ Cotinine and D₄ Nicotine are available as deuterated internal standards.

SOURCE- UCT Internal Publication

**OPIATES IN URINE-OXIME TMS PROCEDURE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)
1 of 2

1. PREPARE SAMPLE

ACID HYDROLYSIS OF GLUCURONIDES:

To 2 mL of urine add internal standard(s)* and 400 µl concentrated HCl.

Add 200 µl 10% Hydroxylamine Solution.

Mix/vortex.

Heat to 90° C for 20 min in a heating block or an autoclave for 15 minutes on a liquid cycle

Cool before proceeding.

Centrifuge for 10 minutes at 2000 rpm and discard pellet.

Add 500 µl 50% Ammonium Hydroxide Mix/vortex.

Adjust sample pH 6.0-7.0 by drop wise addition with 50% Ammonium Hydroxide

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL methanol; aspirate.

1 x 3 mL DI H₂O; aspirate.

1 x 3 mL 0.1 M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.

1 x 3 mL 0.1 M acetate buffer (pH 4.5); aspirate.

1 x 3 mL methanol; aspirate.

Dry column (5 minutes at > 10 inches Hg).

5. ELUTE OPIATES

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2)

Collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C

**OPIATES IN URINE-OXIME TMS PROCEDURE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

7. DERIVATIZE

Add 100 µL ethyl acetate and 50 µL MSTFA .
Overlayer with N₂ and cap. Mix/vortex.
React 20 minutes at 85° C. in a heat block.
Remove from heat source to cool

NOTE: Do not evaporate MSTFA solution.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Quant Ion	Secondary	Tertiary
Codeine TMS	371	356	343
D ₃ -Codeine TMS*	374	359	346
Morphine TMS	236	414	429
D ₃ -Morphine TMS*	239	417	432
Hydrocodone Oxime TMS	386	371	329
Hydromorphone Oxime TMS	355	429	444
Oxycodone-Oxime TMS	474	459	401
Oxморphone-Oxime TMS	459	287	532

*Suggested internal standards for GC/MS: D₃-Codeine, D₃-Morphine

**OPIATES IN HUMAN URINE- PROPYL DERIVATIVES
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1 of 2

- 1. PREPARE SAMPLE-ENZYMATIC HYDROLYSIS OF GLUCURONIDES:**
To 2 mL of urine, add internal standard(s)*, and 1 mL of β -Glucuronidase solution. (β -Glucuronidase solution contains 5,000 F units/mL *Patella Vulgata* in 0.1 M acetate buffer, pH 5.0). Hydrolyze for 3 hours at 60° C. Cool before proceeding. Centrifuge for 10 minutes at 2000 rpm and discard pellet. Adjust sample pH to 6.0 \pm 0.5 with 1.0 N NaOH.
- 2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN**
1 x 2 mL methanol; aspirate.
1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 0.1M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.
- 3. APPLY SAMPLE**
Load at 1 to 2 mL/minute.
- 4. WASH COLUMN**
1 x 3 mL DI H₂O; aspirate.
1 x 3 mL 0.1 M acetate buffer (pH 4.5); aspirate.
1 x 3 mL methanol; aspirate.
Dry column (5 minutes at > 10 inches Hg).
- 5. ELUTE OPIATES**
1 x 2 mL ethyl acetate/isopropanol/ammonium hydroxide (84/12/4)
- 6. DRY ELUANT**
Evaporate to dryness at < 40° C
- 7. DERIVATIZE**
Add 200 μ L of a 1:1 solution of propionic anhydride: pyridine
Make this solution fresh daily.
Mix/vortex.
React for 60 minutes at 40° C in a heater block.
Remove from heat source to cool.
Evaporate to dryness at < 40° C
Reconstitute the residue with 50 μ L of an ethyl acetate / methanol (70:30).

**OPIATES IN HUMAN URINE- PROPYL DERIVATIVES
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound (Propyl)	Primary Ion**	Secondary Ion	Tertiary Ion
Hydrocodone***	299	242	214
Codeine	355	282	229
d ₃ -Codeine*	358	285	232
Oxycodone	371	314	298
Hydromorphone	285	341	228
6-Monoacetylmorphine	327	268	383
Oxymorphone	357	300	413
Morphine	341	268	397
d ₃ -Morphine*	344	271	400

* Suggested internal standard for GC/MS: d₃-Codeine and d₃-Morphine

** Quantitation ion

*** Hydrocodone does not derivatize under these conditions.

SOURCE- UCT Internal Publication working with the Philadelphia Medical Examiner's Office

**FREE (UNBOUND) OPIATES
IN SERUM, PLASMA OR WHOLE BLOOD
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 1 mL of sample (serum, plasma or whole blood) add internal standard(s)* and 4 mL of DI H₂O. Mix/vortex and let stand 5 minutes.
Centrifuge for 10 minutes at 2000 rpm and discard pellet.
Add 2 mL of 0.1 M phosphate buffer (pH 6.0). Mix/vortex.
Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₂OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 0.1M acetate buffer (pH 4.5); aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at >10 inches Hg).

5. ELUTE OPIATES

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 to 2 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C.

7. DERIVATIZE

Add 50 µL ethyl acetate and 50 µL BSTFA (with 1% TMCS).
Overlay with N₂ and cap. Mix/vortex.
React 20 minutes at 70° C. Remove from heat source to cool.
NOTE: Do not evaporate BSTFA solution.

8. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound (TMS)	Primary Ion**	Secondary Ion	Tertiary Ion
Codeine:	371	234	343
D ₃ -Codeine*	374	237	346
Morphine	429	287	324
D ₃ -Morphine*	432	290	327

* Suggested internal standards for GC/MS: D₃-Codeine, D₃-Morphine

** Quantitation ion

PHENCYCLIDINE IN URINE FOR GC OR GC/MS CONFIRMATIONS USING 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine add internal standard(s)* and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix / vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL methanol; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate
1 x 3 mL methanol; aspirate.
Dry column (5 minutes at > 10 inch Hg).

5. ELUTE PHENCYCLIDINE

1 x 3 mL Methylene Chloride/Isopropanol/Ammonium Hydroxide (78/20/2)
NOTE: Prepare elution solvent daily. Add IPA/NO₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C
Remove immediately upon completion.
Reconstitute with 100 µL ethyl acetate.

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

COMPOUND	Primary Ion**	Secondary	Tertiary
Phencyclidine	200	91	242
D ₅ - Phencyclidine*	205	96	247

*Suggested internal standard for GC/MS: D₅-Phencyclidine

**Quantitation Ion

REFERENCE- UCT Internal Publication

PROPOXYPHENE IN URINE FOR GC OR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 2 mL of urine add internal standard(s)* and 1 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex. Sample pH should be 6.0 ± 0.5.
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 2 mL 0.1 M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE PROPOXYPHENE

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 mL/minute.
NOTE: Prepare elution solvent daily. Add IPA/NO₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. CONCENTRATE

Evaporate to dryness at < 40° C.
Reconstitute with 100 µL ethyl acetate.

7. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

Compound	Primary**	Secondary	Tertiary	OTHER
Propoxyphene	58	115	208	250,265
Propoxyphene D ₅ *	63	120	213	255,270

* Internal Standard

** Quantitation Ion

NOTE: To improve the analysis for Norpropoxyphene, the primary metabolite of Dextropropoxyphene, add 1 drop of 35% sodium hydroxide solution to the urine sample and then after mixing bring the pH to 6 for SPE extraction. This step converts the Norpropoxyphene to Norpropoxyphene amide, a more stable compound.

For more information see the following reference:

Amalfitano G, Bessard J, Vincent F, Esseric H and Bessard G
Gas Chromatographic Quantitation of Dextropropoxyphene and Norpropoxyphene in Urine after Solid Phase Extraction
Journal Analytical Toxicology 20:547-554 (1996)

**PSILOCIN IN URINE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1 of 2

1. PREPARE SAMPLE

To 5 mL of urine add internal standard and 2 mL of 0.1 M phosphate buffer (pH 6.0).
Mix/vortex.
Add 12,500 to 25,000 units of β -Glucuronidase, Mix/Vortex
Place the sample into a water bath at 45°C for 90 minutes.
Remove from the bath and allow to cool
Centrifuge at 3,000 rpm for 10 min
Use the clear filtrate (discard the plug) for SPE

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 2 mL DI H₂O; aspirate.
1 x 2 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 2 mL 20% Acetonitrile in water; aspirate
1 x 1 mL 0.1M Acetic Acid; aspirate
Dry column (3 minutes at > 10 inches Hg).
1 x 2 mL Hexane; aspirate.
1 x 3 mL Hexane/ Ethyl Acetate (50:50)
1 x 3 mL Methanol
Dry Column (3 min at > 10 inches Hg)

5. ELUTE PSILOCIN

1 x 3 mL Dichloromethane/ Isopropanol/ Ammonium Hydroxide (78:20:2)
Collect Eluant at 1mL /min
NOTE: Prepare elution solvent daily

6. DRY ELUATE

Evaporate to dryness at < 35°C

7. DERIVATIZE

Reconstitute with 50 μ L MSTFA
Cap the sample tube and place it into a heater block at 70°C for 20 minutes.

**PSILOPIN IN URINE
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

8. QUANTITATE

Inject 1 to 2 µL onto chromatograph.
Monitor the following ions (Mass Selective Detection):

Compound	Primary*	Secondary	Tertiary
PSILOPIN – TMS	290	348	73 (291)

* Quantitation Ion

GC CONDITIONS:

HP Model 5890 GC with a 5970 MSD
COLUMN = DB5 (25M x 0.32mmID x 0.17µ Film Thickness)
CARRIER GAS -Helium (5psi head pressure)

INJECTION Size = 1
| SPLITLESS MODE

Injection Temperature = 275°C
Detector Temperature = 300°C
TEMPERATURE PROGRAM:
70° C hold 1 min then ramp to 240° C at 20°C/min hold for 2 minutes

SOURCE - The Detection of Psilocin in Human Urine
Grieshaber A, Moore K, Levine B and Smith M
Presented at the TRI-SERVICES Meeting Nov 1999

**SERTRALINE AND DESMETHYLSERTRALINE
IN SERUM, PLASMA OR WHOLE BLOOD
FOR HPLC ANALYSIS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE

To 1 mL of sample (serum, plasma or whole blood) add internal standard, 4 mL DI H₂O and 2 mL of 0.1 M phosphate buffer (pH 6.0).

Mix/vortex. Centrifuge for 10 minutes at 2000 rpm and discard pellet

Sample pH should be 6.0 ± 0.5.

Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.

1 x 3 mL DI H₂O; aspirate.

1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.

1 x 1 mL 0.1 M acetic acid; aspirate.

1 x 3 mL CH₃OH; aspirate.

Dry column (5 minutes at > 10 inches Hg).

5. ELUTE

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2); collect eluate at 1 mL/minute.

NOTE: Prepare elution solvent fresh daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C

7. QUANTITATE

Reconstitute with 200 µL acetonitrile/DI H₂O (1/3). Mix/vortex vigorously for 30 seconds.

Inject 100 µL onto chromatograph at wavelength 235 nm.

Mobile phase = 0.25M potassium phosphate (pH 2.7) containing 30% CH₃CN

Flow rate = 2 mL/minute.

HPLC SYSTEM:

Isocratic HPLC using a Pump thru a C8 HPLC Column

(LC-8 or equivalent HPLC Column) 15cm X 4.6mm ID

Coupled to a UV detector set at 235 nm.

THC AND CARBOXY THC IN WHOLE BLOOD FOR GC/MS CONFIRMATIONS USING: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSTHC020 without Tips or ZCTHC020 with CLEAN-THRU® Tips)

1 of 2

1. PREPARE SAMPLE

To 1 mL of whole blood sample add internal standard(s)* into a silanized screw top test tube.
Mix/vortex. Let stand up 1 hour.
Add 2 mL of acetonitrile, Vortex for 30 seconds and allow to sit for 3 minutes.
Centrifuge for 10 minutes at maximum rpm.
Decant into another silanized test tube and add 4 mL of 0.1 M acetate buffer (pH 4.5) to supernatant.
Mix/vortex, centrifuge at maximum rpm for 5 minutes and decant the supernatant into another silanized test tube. This step removes any blood fragments or foam.

2. PRECONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL hexane/ethyl acetate (75/25); aspirate.

3. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.

1 x 3 mL DI H₂O; aspirate.

1 x 1 mL 0.1 M HCl; aspirate.

NOTE: Use gravity flow or minimal vacuum.

4. APPLY SAMPLE

Load at 1 mL/minute.

NOTE: Use gravity flow or minimal vacuum.

5. WASH COLUMN

1 x 2 mL DI H₂O; aspirate.

1 x 2 mL 0.1 M HCl/acetonitrile (70/30); aspirate.

Dry columns (5 minutes at > 10 inches Hg).

1 x 200 µL hexane.

NOTE: Use gravity flow or minimal vacuum.

6. ELUTE THC AND CARBOXY THC

1 x 3 mL hexane/ethyl acetate (75/25).

NOTE: Use gravity flow or minimal vacuum.

7. DRY ELUATE

Evaporate slowly to dryness at < 40° C.

8. DERIVATIZE

Add 50 µL BSTFA (with 1% TMCS) and 50 µL of ethyl acetate

Overlayer with Nitrogen gas and cap.

Mix/vortex.

React 30 minutes at 70° C.

Remove from heat source to cool.

NOTE: Do not evaporate BSTFA solution.

**THC AND CARBOXY THC IN WHOLE BLOOD
FOR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

9. QUANTITATE

Inject 1 to 2 µL onto gas chromatograph.
For MSD monitor the following ions:

ANALYTE (TMS)	Primary Ion**	Secondary	Tertiary
Carboxy- Δ^9 THC	371	473	488
D ₃ Carboxy- Δ^9 THC*	374	476	491
Tetrahydrocannabinol	371	303	386
D ₃ Tetrahydrocannabinol *	374	306	389

*Suggested internal standards for GC/MS: D₃THC and D₃Carboxy- Δ^9 THC

**Quantitation ion

Reference- UCT Internal Publication

**THERAPEUTIC AND ABUSED DRUGS IN URINE
FOR ACID/NEUTRAL AND BASIC DRUGS
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN[®] EXTRACTION COLUMN**

**(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU[®] Tips)
1 of 3**

1. PREPARE SAMPLE

Urine

To 2 mL of urine add internal standard(s) and 1 mL of 0.1 M phosphate buffer pH 6.0.
Mix/vortex. Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

Serum, Plasma or Whole Blood

To 1 mL of sample add internal standard(s) and 4 mL DI H₂O (5.5 pH 5.7).
Mix/vortex and let stand 5 minutes.
Centrifuge for 10 minutes at 2000 rpm and discard pellet.
Add 2 mL 0.1 M phosphate buffer (pH 6.0).
Mix/vortex.
Sample pH should be 6.0 ± 0.5 .
Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN[®] EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
Dry column (5 minutes at > 10 inches Hg).
1 x 2 mL hexane; aspirate.

5. ELUTE ACIDIC AND NEUTRAL DRUGS (FRACTION 1)

1 x 3 mL hexane/ethyl acetate (50/50); collect eluate at < 2 mL/minute.

6. DRY ELUATE

Evaporate to dryness at < 40° C
Reconstitute with 100 µl ethyl acetate.

7. QUANTITATE ACIDIC AND NEUTRAL DRUGS

Inject 1 to 2 µl onto gas chromatograph.

8. WASH COLUMN

1 x 3 mL methanol; aspirate.
Dry column (5 minutes at > 10 inches Hg).

**THERAPEUTIC AND ABUSED DRUGS IN URINE
FOR ACID/NEUTRAL BASIC DRUGS
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 3

9. ELUTE BASIC DRUGS (FRACTION 2)

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2);

Collect eluate at 1 to 2 mL/minute.

NOTE: Prepare elution solvent fresh daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

10. DRY ELUATE

Evaporate to dryness at < 40° C using a TurboVap® or equivalent evaporator.

Take care not to overheat or over-evaporate.

Certain compounds are heat labile, such as the amphetamines and phencyclidine.

Reconstitute with 100 µL methanol or ethyl acetate.

11. QUANTITATE Basic Drugs

Inject 1 to 2 µL onto gas chromatograph.

NOTES:

- (1) Fraction 1 (Acid Neutrals) and Fraction 2 (Bases) can be combined together.
- (2) A keeper solvent such as DMF can be used to prevent the volatilization of amphetamines and phencyclidine. Use 30-50 µL of high purity DMF in the sample (Fraction 2) before evaporation.
- (3) A 1% HCl in MeOH solution has been used to prevent volatilization by the formation of the hydrochloric salt of the drugs. Evaporate fraction 2 to approximately 100 µL, then add 1 drop of the solution. Continue to evaporate to dryness.

SOURCE: UCT Internal Publication

CLEAN SCREEN® DAU Forensic Applications

3 of 3

Data Provided By:

City of Philadelphia, Department of Public Health
Office of the Medical Examiner
321 University Avenue
Philadelphia, Pennsylvania 19104

Contact: Frank Caputo, Analytical Chemist II
(215) 823-7464

The following are some of the many compounds that have been extracted from forensic samples with the CLEAN SCREEN® DAU bonded silica extraction cartridge (**part # CSDAU303**):

I. ACIDIC/NEUTRAL DRUG FRACTION (A)

Acetaminophen	Clonazepam	Nordiazepam
Barbiturates	Cotinine	Phenytoin
Benzoic acid	Diazepam	Primidone
Caffeine	Glutethimide and metabolite	Salicylic acid
Carbamazepine	Ibuprofen	Theophylline
Carisoprodol	Meprobamate	Thiopental
Chlorpropamide	Methyl salicylate	

II. BASIC DRUG FRACTION (B)

Amantadine	Dihydrocodeine	Methylphenidate
Amitriptyline and metabolite	Diethylhydramine	Methyprylon and metabolite
Amphetamine	Doxepin and metabolite	Morphine
Benzocaine	Ephedrine	Nicotine
Benzoylcegonine	Fluoxetine	Oxycodone
Benztropine	Imipramine and metabolite	Pentazocine
Bromodiphenhydramine	Ketamine	Phencyclidine
Chlordiazepoxide	Lidapine	Phenethylamine
Chloroquine	Loxapine	Phentermine
Chlorpheniramine	Meperidine	Phenylpropanolamine
Chlorpromazine	Methadone and metabolite	Procaine
Cocaine and metabolite	Methamphetamine	Propoxyphene and metabolite
Codeine	Methyl <i>p</i> -aminobenzoate	Propylparaben
Cresol	Methyl benzoate	Tranlycypromine
Dextromethorphan	Methyl ecgonine	Trifluoperazine
Dextrorphan	Methylparaben	Trimipramine
		Thioridazine
		Trazodone

United Chemical Technologies, Inc.
2731 Bartram Road
Bristol, PA 19007
USA

(215) 781-9255
Toll Free: (800) 541-0559
Fax: (215) 785-1226

**TRICYCLIC ANTIDEPRESSANTS IN SERUM AND PLASMA
FOR HPLC USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1 of 2

1. PREPARE SAMPLE

To 1 mL of serum or plasma add internal standard and 2 mL of 0.1 M phosphate buffer (pH 6.0). Mix/vortex. Centrifuge for 10 minutes at 2000 rpm and discard pellet. Sample pH should be 6.0 ± 0.5 . Adjust pH accordingly with 0.1 M monobasic or dibasic sodium phosphate.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN

1 x 3 mL CH₃OH; aspirate.
1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE

Load at 1 mL/minute.

4. WASH COLUMN

1 x 3 mL DI H₂O; aspirate.
1 x 1 mL 0.1 M acetic acid; aspirate.
1 x 3 mL CH₃OH; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE

1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78/20/2)
Collect eluate at 1 mL/minute or use gravity flow
NOTE: Prepare elution solvent fresh daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. DRY ELUATE

Evaporate to dryness at < 40° C

7. QUANTITATE

Reconstitute with 200 µl acetonitrile/DI H₂O (1/3).
Mix/vortex vigorously for 30 seconds.
Inject 100 µL onto HPLC

**TRICYCLIC ANTIDEPRESSANTS IN SERUM AND PLASMA
FOR GC OR GC/MS CONFIRMATIONS USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**

2 of 2

HPLC CONDITIONS

HPLC COLUMN – Propylcyano, Endcapped 4.6mm x 150mm, 5µ particle size
 COLUMN TEMPERATURE = 30° C
 MOBILE PHASE- Acetonitrile/ Buffer/ Methanol (60/25/15)
 Buffer= 0.01M K₂HPO₄ adjusted to pH 7.0 with H₃PO₄
 FLOW RATE = 1.75 mLs/min

ANALYTES AND EXTRACTION EFFICIENCY

COMPOUND	Retention Time (min)	% Recovery	%RSD
Trimipramine ISTD*	2.564	100.0%	5.53%
Doxepin	3.048	96.5%	8.04%
Amitriptyline	3.433	98.9%	5.64%
Imipramine	3.865	97.2%	6.09%
Nortriptyline	5.349	88.9%	9.49%
Nordoxepin	5.788	85.0%	5.29%
Desipramine	6.067	85.3%	5.04%
Protriptyline ISTD*	6.476	86.3%	5.39%

* Internal Standards

HINTS:

- (1) Silica Based HPLC columns are sensitive to pH.
 To prevent dissolution of the packing especially at the head of the column, it is best to place a silica column before the injector. This will saturate the mobile phase with silica.
- (2) Secondary Amines bind to glass and polyethylene.
 It is recommended to silylate all surfaces that come in contact with the sample.
 Immersion into 5% DMCS in toluene or vapor deposition will deactivate the surface by silylation.
- (3) To ensure the proper strength of Elution solvent.
 Measure the apparent pH of the elution solvent.
 It should be pH 10 or better. Add 1-2% of Ammonium Hydroxide and check again.

**4'4-METHYLENEDIANALINE
IN SERUM
(ENVIRO-CLEAN® C18 BONDED PHASE)**

(Part # CEC18111)

- 1. SAMPLE PREPARATION**
 - a. none
- 2. CONDITION ENVIRO-CLEAN® EXTRACTION COLUMN**
 - a. 3 mL Methanol
 - b. 3 mL Distilled Water
- 3. LOAD 1 mL SERUM SAMPLE**
- 4. COLUMN WASH**
 - a. 1 mL Distilled Water
- 6. ELUTION**
 - a. 0.25 mL Methanol containing 1M ammonium hydroxide
- 7. INJECTION / ANALYSIS**
 - a. Inject 10 µL onto HPLC system

4'4 Methylenedianaline
Extracted 100 µg/mL sample
Retention time = 3.162 min
Mobile Phase: Methano/H₂O (50:50)
Selectra C18 column
Flow rate = 1.2 mL/m
Inj. volume = 10 µL
wavelength = 254 nm

**CHLOROPHENOXY ACID HERBICIDES
IN WATER
(ENVIRO-CLEAN® C18 BONDED PHASE)**

(Part # CEC181M6)

Chlorophenoxy Acid Herbicides Extracted
2, 4-D Acid 2,4,5-trichloro phenoxy propionic acid (Silvex) Dicamba Dinitro-sec-butyl phenol

- 1. SAMPLE PREPARATION**
 - a. Adjust pH of 1 liter water sample to pH 1.0 with hydrochloric acid
- 2. CONDITION ENVIRO-CLEAN® EXTRACTION COLUMN**
 - a. 10 mL Hexane / Acetone (50:50)
 - b. 10 mL Acidified Methanol (5% HCl in MEOH)
 - c. 10 mL Distilled water
- 3. LOAD 1 LITER OF PH ADJUSTED WATER SAMPLE (RATE 8-10 ML/MIN)**
- 4. COLUMN WASH**
 - a. 10 mL Distilled water adjusted to pH 1.0 with HCl
- 5. DRY COLUMN FOR 15-30 MINUTES UNDER MAXIMUM VACUUM PRESSURE**
- 6. ELUTION**
 - a. 10 mL Hexane / Acetone (50:50)
- 7. CONCENTRATION / EVAPORATION**
 - a. Add 500 µL of a keeper solvent (Methanol, DMF, other). Evaporate to 500 µL under a nitrogen stream at room temperature
- 8. INJECTION / ANALYSIS**
 - a. Reconstitute with 100 µL TCTEF Inject at 1-2 µL onto GC

**PESTICIDES
IN WATER
(ENVIRO-CLEAN® C18 BONDED PHASE)**

(Part # CEC181M6)

Pesticides Extracted	
α- hexachlorocycchexane	4,4'-DDE
Lindane	Dieldrin
β-hexachlorocycchexane	Endrin
Heptechlor	4,4'-DDD
δ-hexachlorocycchexane	Endosulfan II
Aldrin	4,4'-DDT
Heptachlor	Endrin Aldehyde
Endosulfan 1	Endosulfan Sulfate

1. **SAMPLE PREPARATION**
 - a. None

2. **CONDITION ENVIRO-CLEAN® EXTRACTION COLUMN**
 - a. 10 mL Hexane / Acetone (50:50)
 - b. 10 mL Methanol
 - c. 10 mL Distilled Water

3. **LOAD 1 LITER SAMPLE (RATE 8-10 ML/MIN)**

4. **COLUMN WASH**
 - a. 20 mL Distilled Water

5. **DRY COLUMN FOR 15-30 MINUTES UNDER MAXIMUM VACUUM PRESSURE**

6. **ELUTION**
 - a. 10 mL Hexane / Acetone (50:50)

7. **CONCENTRATION / EVAPORATION**
 - a. Add 500 μL of a keeper solvent (methanol, DMF, other).
 - b. Evaporate to 500 μL under a nitrogen stream

**Polychlorinated Biphenyls
In Pond Water
(ENVIRO-CLEAN® C18 Bonded Phase)
(Part # CEC182M6)**

Polychlorinated Biphenyls Extracted	
Aroclor 1026	Aroclor 1221
Aroclor 1232	Aroclor 1242
Aroclor 1248	Aroclor 1254
Aroclor 1260	

- 1. SAMPLE PREPARATION**
 - a. Filter water sample through a 0.5 µm filter.
 - b. Add 2 mL methanol to 200 mL of filtered water sample.
 - c. Mix and degas sample for 2 minutes

- 2. CONDITION ENVIRO-CLEAN® EXTRACTION COLUMN**
 - a. 15-20 mL Hexane Dry for 5 minutes
 - b. 15-20 mL Methanol
 - c. 10 mL Distilled Water

- 3. LOAD 200 ML WATER SAMPLE (RATE 8-10 ML/MIN)**

- 4. COLUMN WASH**
 - a. 20 mL Distilled Water

- 5. DRY COLUMN FOR 20-30 MINUTES UNDER MAXIMUM VACUUM PRESSURE**

- 6. ELUTION**
 - a. 20 mL Hexane

- 7. CONCENTRATION / EVAPORATION**
 - a. Evaporate to dryness under a nitrogen stream at room temperature

- 6. INJECTION / ANALYSIS**
 - a. Reconstitute with 100 µL methanol inject at 1-2 µL aliquot onto GC

**POLYNUCLEAR AROMATIC HYDROCARBONS
IN POND WATER
(ENVIRO-CLEAN® C18 BONDED PHASE)**

(Part # CUNAX153)

Polynuclear Aromatic Hydrocarbons Extracted	
Naphthalene	Chrysene
Fluorene	B(e)pyrene
Acenaphthene	B(b)fluoranthene
Phenanthrene	B(k)fluoranthene
Anthracene	B(a)pyrene
Fluoranthene	D(a,h)anthracene
Pyrene	B(g,hi)perylene
B(a)anthracene	Indeno(1,2,3,-cd)pyrene

- 1. SAMPLE PREPARATION**
 - a. Filter water through a 0.5 µm filter.
 - b. Add 2 mL of methanol to 200 ml of filtered water sample.
 - c. Mix and degas sample for 2 minutes.

- 2. PREPARE ENVIRO-CLEAN® EXTRACTION COLUMN**
 - a. 15-20 mL Methylene chloride/trichlorotrifluoroethylene (TCTFE)
 - b. 15-20 mL TCTFE. Dry for 5 minutes
 - c. 15-20 mL Methanol
 - d. 20 mL Distilled water

- 3. LOAD 200 mL WATER SAMPLE (RATE 8-10 ML/MIN)**

- 4. COLUMN WASH**
 - a. 20 ml Distilled Water

- 5. DRY COLUMN FOR 15-30 MINUTES UNDER MAXIMUM VACUUM PRESSURE**

- 6. ELUTION**
 - a. 20 mL of TCTFE

- 7. CONCENTRATION / EVAPORATION**
 - a. Evaporate to dryness under a nitrogen stream at room temperature

- 8. INJECTION / ANALYSIS**
 - a. Reconstitute with 100 µL TCTEF Inject at 1-2 µL aliquot onto GC

ABUSED DRUGS IN CANINE OR EQUINE URINE USING: 500 mg XtrackT[®] EXTRACTION COLUMN

(Part # XRDAH515)

1. PREPARE SAMPLE-ENZYMATIC HYDROLYSIS OF GLUCURONIDES

To 5 mL of urine add internal standard(s) and 2 mL of β -Glucuronidase 5,000 F units/mL *Patella Vulgata* in 100 mM Acetate Buffer (pH 5.0).

Mix/vortex. Hydrolyze at 65°C for 3 hours.

Cool before proceeding.

BASE HYDROLYSIS OF GLUCURONIDES

To 2 mL of urine add internal standard(s) and 100 μ L of 10N NaOH

Mix/vortex. Hydrolyze at 60°C for 20 minutes.

Cool before proceeding.

COMBINE HYDROLYSATES

Combine both hydrolysis products with 5 mLs of 0.1M pH 6.0 Phosphate Buffer

Adjust sample pH = 6.0 \pm 0.5 with 0.5M Phosphoric acid

2. CONDITION XtrackT[®] EXTRACTION COLUMN

1 x 5 mL CH₃OH; aspirate.

1 x 5 mL DI H₂O; aspirate.

1 x 3 mL 0.1 M phosphate buffer (pH 6.0); aspirate.

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE Load at 1 to 2 mL/minute.

4. WASH COLUMN

1 x 3 mL 0.1M pH 6.0 Phosphate buffer; aspirate.

1 x 2 mL 1.0 M acetic acid; aspirate.

Dry column (5 minutes at > 10 inches Hg).

1 x 2 mL hexane; aspirate.

5. ELUTE ACIDIC AND NEUTRAL DRUGS

1 x 4 mL methylene chloride; collect eluate at < 5 mL/minute.

6. ELUTE STEROIDS 2 x 4 mL ethyl acetate; collect eluate at < 5 mL/minute.

7. WASH COLUMN 1 x 5 mL methanol; aspirate

8. ELUTE BASIC DRUGS

1 x 5 mL methylene chloride/ isopropanol/ ammonium hydroxide (78/20/2)

NOTE: Prepare elution solvent fresh daily

9. DRY ELUATE

Evaporate to dryness at < 40° C

Reconstitute with 100 μ L ethyl acetate.

10. QUANTITATE Spot onto TLC plate or inject 1 to 2 μ L onto chromatograph

EXTRACTION OF TEAR GAS

**Chloroacetophenone (CS), o-Chlorobenzylidenemalononitrile (CN),
and trans-8-methyl-N-vanillyl-6-nonenamide (OC)**

From Cloth for GC/MS Analysis Using: 200 mg CLEAN SCREEN® EXTRACTION COLUMN

(Part # ZSDAU020 without Tips or ZCDAU020 with CLEAN-THRU® Tips)

1. PREPARE SAMPLE:

If suspected tear gas is on clothing cut out a portion of the sprayed area and a “negative” control sample. Extract each into hexane.
For canisters of suspected tear gas, spray onto a Kimwipe® and extract the sprayed area and a negative control into hexane.

2. CONDITION CLEAN SCREEN® EXTRACTION COLUMN:

2 x 3 mL CH₃OH; aspirate.
2 x 3 mL DI H₂O; aspirate.
1 x 2 mL 100 mM phosphate buffer (pH 6.0); aspirate.
NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. APPLY SAMPLE:

Load at 1 ml/minute.

4. WASH COLUMN:

1 x 3 mL DI H₂O; aspirate.
1 x 3 mL Hexane; aspirate.
Dry column (5 minutes at > 10 inches Hg).

5. ELUTE ANALYTE:

1 x 1 mL CH₃OH.

6. DRY ELUATE:

Evaporate to dryness at < 40° C.

7. RECONSTITUTE

Add 200 µL CH₃OH. Mix/vortex. Transfer to GC/MS vial and cap.

8. QUANTITATE

Inject 1-2 µL sample onto GC/MS.

GC/MS Conditions:

Column: HP Ultra 1, Crosslink Methyl Silicone
12 mm x 0.2 x 0.33 µm film thickness

GC Oven: Initial Temp. = 100° C
Initial Time = 3.00 min.

Ramp = 17C/min.
Final Temp. = 305° C
Final Time = 3.00 min.
Injection Port Temp. = 250° C
Transfer line Temp. = 280° C

SCAN Acquisition = 41 amu to 400 amu: Start time = 2.00 mins

Retention times:	Compound	CN	CS	OC
	RT (min.)	@4.9	@7.4	@13.4



UNITED CHEMICAL TECHNOLOGIES, INC.

**Return to
SPE Products
Catalog**

**Return to
SPE Applications
Manual**

**Return to
Research Papers**

Get the **Scoop on the Latest
Solid Phase Extraction Methods!**

Call Your UCT Technical Specialist
Toll Free at 800-541-0559

or

Visit www.unitedchem.com

